

Product datasheet for **RC222564L4V**

Gastrokine 1 (GKN1) (NM_019617) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gastrokine 1 (GKN1) (NM_019617) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Gastrokine 1
Synonyms:	AMP18; BRICD1; CA11; FOV; foveolin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019617
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222564).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019617.3 , NP_062563.3
RefSeq Size:	820 bp
RefSeq ORF:	558 bp
Locus ID:	56287
UniProt ID:	Q9NS71
Cytogenetics:	2p13.3
Protein Families:	Secreted Protein, Transmembrane
MW:	21.8 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is found to be down-regulated in human gastric cancer tissue as compared to normal gastric mucosa. [provided by RefSeq, Jul 2008]