

Product datasheet for **RC222571L4V**

ZNF12 (NM_006956) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF12 (NM_006956) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF12
Synonyms:	GIOT-3; HZF11; KOX3; ZNF325
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006956
ORF Size:	1977 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222571).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006956.2
RefSeq Size:	4960 bp
RefSeq ORF:	1980 bp
Locus ID:	7559
UniProt ID:	P17014
Cytogenetics:	7p22.1
Protein Families:	Transcription Factors
MW:	77 kDa


[View online »](#)

Gene Summary:

This gene is a member of the krueppel C2H2-type zinc-finger protein family and encodes a protein with eight C2H2-type zinc fingers and a KRAB domain. This nuclear protein is involved in developmental control of gene expression. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]