

Product datasheet for **RC222572L3V**

ANKS1B (NM_020140) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ANKS1B (NM_020140) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ANKS1B
Synonyms:	AIDA; AIDA-1; ANKS2; cajalin-2; EB-1; EB1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020140
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222572).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020140.2 , NP_064525.1
RefSeq Size:	1714 bp
RefSeq ORF:	1245 bp
Locus ID:	56899
UniProt ID:	Q7Z6G8
Cytogenetics:	12q23.1
Domains:	SAM, PID
MW:	46.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a multi-domain protein that is predominantly expressed in brain and testis. This protein interacts with amyloid beta protein precursor (AbetaPP) and may have a role in normal brain development, and in the pathogenesis of Alzheimer's disease. Expression of this gene has been shown to be elevated in patients with pre-B cell acute lymphocytic leukemia associated with t(1;19) translocation. Alternatively spliced transcript variants encoding different isoforms (some with different subcellular localization, PMID:15004329) have been described for this gene. [provided by RefSeq, Aug 2011]