

Product datasheet for **RC222575L4V**

WIZ (NM_021241) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WIZ (NM_021241) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WIZ
Synonyms:	ZNF803
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021241
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222575).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021241.2
RefSeq Size:	4185 bp
RefSeq ORF:	2385 bp
Locus ID:	58525
UniProt ID:	O95785
Cytogenetics:	19p13.12
MW:	85.2 kDa



[View online »](#)

Gene Summary:

May link EHMT1 and EHMT2 histone methyltransferases to the CTBP corepressor machinery.
May be involved in EHMT1-EHMT2 heterodimer formation and stabilization (By similarity).
[UniProtKB/Swiss-Prot Function]