

Product datasheet for **RC222586A1V**

Human UBE2D3 (NM_181888) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human UBE2D3 (NM_181888) AAV Particle
Tag: Myc-DDK
Symbol: UBE2D3
Synonyms: E2(17)KB3; UBC4/5; UBCH5C
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC222586 representing NM_181888
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGCTGAAACGGATTAATAAGGAACCTAGTGATTTGGCCCGTGACCCTCCAGCACAAATGTTCTGCAG
GTCCAGTTGGGGATGATATGTTTCATTGGCAAGCCACAATTATGGGACCTAATGCAGCCCATATCAAGG
CGGTGTATTCTTTTTGACAATTCATTTTCCTACAGACTACCCCTCAAACCACCTAAGGTTGCATTTACA
ACAAGAATTTATCATCCAAATATTAACAGTAATGGCAGCATTTGTCTCGATATTCTAAGATCACAGTGGT
CGCTGCTTTAAACAATTTCTAAAGTTCTTTTATCCATTTGTTCACTGCTATGTGATCCAAACCCAGATGA
CCCCCTAGTGCCAGAGATTGCACGGATCTATAAAACAGACAGAGATAAGTACAACAGAATATCTCGGAA
TGGACTCAGAAGTATGCCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222586 representing NM_181888
Red=Cloning site Green=Tags(s)

MALKRINKELSDLARDPPAQCSAGVPVGDMMFHWQATIMGPNDSPYQGGVFFLTIHFPTDYPFKPPKVAFT
TRITYHPNINSNGSICLDILRSQWSPAL TISKVLLSICSLLCDPNPDDPLVPEIARIYKTRDKYNRISRE
WTQKYAM

TRTRP**LEQKLI**SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online >](#)

ACCN:	NM_181888
ORF Size:	443 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_181888.1</u>
RefSeq Size:	3833 bp
RefSeq ORF:	444 bp
Locus ID:	7323
UniProt ID:	<u>P61077</u>
Cytogenetics:	4q24
MW:	16.7 kDa