

Product datasheet for **RC222588A1V**

Human OR2T3 (NM_001005495) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR2T3 (NM_001005495) AAV Particle
Tag:	Myc-DDK
Symbol:	OR2T3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC222588 representing NM_001005495 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGCTCAGGGAATCAGACTTCTCAGAATCAAACAGCAAGCACTGATTCACCCTCACGGGACTCTTTG
CTGAGAGCAAGCATGCTGCCCTCCTCTACACCGTGACCTTCTTTCTTGATGGCCCTCACTGGGAA
TGCCCTCCTCATCCTCCTCATCCACTCAGAGCCCCGCTCCACACCCCATGTACTTCTTCATCAGCCAG
CTCGCGCTCATGGATCTCATGTACCTATGCGTGACTGTGCCAAGATGCTTGTGGGCCAGGCTCACTGGAG
ATGATACCATTTCCCGTCAGGCTGTGGGATCCAGATGTTCTTCTACCTGACCCTGGCTGGAGCTGAGGT
TTTCCTCCTGGCTGCCATGGCCTATGACCGATATGCTGCTGTTTGCAGACCTCTCCATTACCCACTGCTG
ATGAACCAAGAGGTGTGCCAGCTCCTGGTGTGACGCTGCTGGGTTTTGGGAATGGTTGATGGTTTTGTC
TCACCCCATTAACCATGAGCTTCCCTTTTGCCAGTCTAGGAAAATCCTGAGTTTTTCTGTGAGACTCC
TGCCCTGCTGAAGCTCTCCTGCTCTGACGTCTCCCTCTATAAGACGCTCATGTACCTGTGCTGCATCCTC
ATGCTTCTCGCCCCATCATGGTCATCTCCAGCTCATACACCTCATCTGCATCTCATCCACAGGATGA
ATTCTGCCGCCGCGCCACAGGAAGGCCTTGCCACCTGCTCCTCCACATGATCATAGTGCTGCTGCTCTT
CGGTGCTTCCTTCTACACCTACATGCTCCCGAGTTCCTACCACACAGCTGAGCAGGACATGATGGTGTCT
GCCTTTTACACCATTTCACTCCTGTGCTGAACCCCTCATTTACAGTCTCCGCAACAAAGATGTACCA
GGGCTCTGAGGAGCATGATGCAGTCAAGAATGAACCAAGAAAAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC222588 representing NM_001005495 Red=Cloning site Green=Tags(s) MCSGNQTSQNQTASTDFTLTGLFAESKHAALLYTVTFLLFLMALTGNALLILLIHSEPRLHTPMYFFISQ LALMDLMYLCVTPKMLVGQVTGDDTISPSGCGIQMFFYLTLAGA EVFLLAAMAYDRYA AAVCRPLHYPLL MNQ RVCQLLV SACWVLGMVDG LLLTPITMSFPFCQSRKILSFFCETPALLKLSCSDVSLYKTLMYLCCIL MLLAPIMVISSSYTLILHLIHRMNSAAGHRKALATCSSHMIIVLLLF GASFYTYMLPSSYHTAEQDMMVS AFYTI FTPVLNPLIYSLRNKDVTRALRSMQSRMNQEK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001005495
ORF Size:	954 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001005495.1, NP_001005495.1</u>
RefSeq Size:	957 bp
RefSeq ORF:	957 bp
Locus ID:	343173
UniProt ID:	<u>Q8NH03</u>
Cytogenetics:	1q44
MW:	35.4 kDa