

Product datasheet for **RC222596A1V**

Human PLGLB1 (PLGLB2) (NM_002665) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PLGLB1 (PLGLB2) (NM_002665) AAV Particle
Tag: Myc-DDK
Symbol: PLGLB1
Synonyms: PLGLA; PLGLA1; PLGP1; PLGP2; PRGA
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC222596 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAACATAAGGAAGTGGTTCTTCTACTTCTTTTATTCTGAAATCAGGTCAAGGAGAGCCCCTGGATG
 ACTATGTGAATACCCAGGGCCCTTCACTGTTCACTGTCAGTGTCACTAAGAAGCAGCTGGGGCAGGAAGCAGAGA
 AGAATGTGCAGCAAAATGTGAAGAGGACAAAGAATTCACCTGCAGGGCATTCCAATATCACAGTAAAGAG
 CAACAATGTGTGATAATGGCTGAAAACAGGAAGTCTCCATAATCATTAGGATGAGAGATGCAGTTTTAT
 TTGAAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222596 protein sequence
 Red=Cloning site Green=Tags(s)

MEHKEVLLLLLLFLKSGQGEPLDDYVNTQGPSLFSVTKKQLGAGSREECAAKCEEDKEFTCRAFQYHSKE
 QQCVIMAENRKSSIIIRMDAVLFEK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_002665
ORF Size: 288 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002665.3</u>
RefSeq Size:	3016 bp
RefSeq ORF:	291 bp
Locus ID:	5342
UniProt ID:	<u>Q02325</u>
Cytogenetics:	2p11.2
MW:	11 kDa