

Product datasheet for **RC222601L3V**

TMEM16B (ANO2) (NM_020373) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM16B (ANO2) (NM_020373) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM16B
Synonyms:	C12orf3; TMEM16B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020373
ORF Size:	2997 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222601).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020373.1 , NP_065106.1
RefSeq Size:	3717 bp
RefSeq ORF:	2996 bp
Locus ID:	57101
Cytogenetics:	12p13.31
Protein Families:	Transmembrane
MW:	113.3 kDa



[View online »](#)

Gene Summary: ANO2 belongs to a family of calcium-activated chloride channels (CaCCs) (reviewed by Hartzell et al., 2009 [PubMed 19015192]).[supplied by OMIM, Jan 2011]