

Product datasheet for **RC222640L3V**

C2CD4B (NM_001007595) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C2CD4B (NM_001007595) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C2CD4B
Synonyms:	FAM148B; NLF2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007595
ORF Size:	1092 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222640).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001007595.1
RefSeq Size:	1329 bp
RefSeq ORF:	1095 bp
Locus ID:	388125
UniProt ID:	A6NLJ0
Cytogenetics:	15q22.2
MW:	38.6 kDa



[View online »](#)

Gene Summary:

May be involved in inflammatory process. May regulate cell architecture and adhesion.
[UniProtKB/Swiss-Prot Function]