

## Product datasheet for **RC222645L4V**

### ZNF498 (ZSCAN25) (NM\_145115) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF498 (ZSCAN25) (NM_145115) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF498
Synonyms:	ZNF498
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145115
ORF Size:	1632 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222645).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_145115.2</a>
RefSeq Size:	4359 bp
RefSeq ORF:	1635 bp
Locus ID:	221785
UniProt ID:	<a href="#">Q6NSZ9</a>
Cytogenetics:	7q22.1
Domains:	zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 61.3 kDa

**Gene Summary:** This gene encodes a protein that bears some similarity to zinc finger proteins, which are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants have been identified, but the full-length nature for most of them has not been determined. [provided by RefSeq, Jul 2008]