

## Product datasheet for **RC222654L1V**

### ASIP (NM\_001672) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | ASIP (NM_001672) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | ASIP   |
| Synonyms:                 | AGSW; AGTI; AGTIL; ASP; SHEP9  |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)  |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001672  |
| ORF Size:                 | 396 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC222654).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001672.2</a> , <a href="#">NP_001663.2</a>  |
| RefSeq Size:              | 584 bp   |
| RefSeq ORF:               | 399 bp   |
| Locus ID:                 | 434  |
| UniProt ID:               | <a href="#">P42127</a>   |
| Cytogenetics:             | 20q11.22   |
| Protein Families:         | Secreted Protein   |
| Protein Pathways:         | Melanogenesis  |



[View online »](#)

**MW:** 14.52 kDa

**Gene Summary:** In mice, the agouti gene encodes a paracrine signaling molecule that causes hair follicle melanocytes to synthesize pheomelanin, a yellow pigment, instead of the black or brown pigment, eumelanin. Pleiotropic effects of constitutive expression of the mouse gene include adult-onset obesity, increased tumor susceptibility, and premature infertility. This gene is highly similar to the mouse gene and encodes a secreted protein that may (1) affect the quality of hair pigmentation, (2) act as a pharmacological antagonist of alpha-melanocyte-stimulating hormone, (3) play a role in neuroendocrine aspects of melanocortin action, and (4) have a functional role in regulating lipid metabolism in adipocytes. [provided by RefSeq, Jul 2008]