

Product datasheet for **RC222658**

TDRKH (NM_001083965) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TDRKH (NM_001083965) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TDRKH
Synonyms:	TDRD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC222658 representing NM_001083965
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCTACTGAACGGACTTCTTGGACAAGCCTGTCCACCATTAGAAAATAGCCCTGGGCCCTGGGATCC
 CAGCCAGTGCAACAGTTGCCTATATCCTATACCGCAGGTATAGGGAAAGCAGAGAAGAGCGGCTGACATT
 TGTGGGGAAGATGACATTGAGATAGAGATGCGGGTCCCCAGGAGGCTGTGAAACTCATCATTGGCCGG
 CAAGGAGCCAATATTAACAGCTGCGGAAACAGACAGGTGCTCGGATTGATGTGGACACAGAGGATGTAG
 GCGATGAGCGAGTGTCTTATCAGTGGTTTTCTGTTCAGGTGTGCAAGGCCAAAGCAGCAATCCATCA
 GATCCTGACAGAGAATACCCAGTGTCTGAGCAGCTTTCAGTCCCCAGAGATCTGTGGCAGAATCATA
 GGGAGAGGCGGAGACAATTCGTTCTATCTGTAAGGCATCTGGAGCCAAAATTACCTGTGACAAAGAAT
 CAGAAGGGACATTACTACTATCAAGACTTATAAAAATCTCAGGAACACAGAAGGAAGTGGCAGCAGCCAA
 GCATTTGATACTGGAGAAAGTTTCAGAAGATGAAGAACTTCGGAAGAGAATTGCTCATTCTGCAGAAACC
 AGGGTCCACGCAACAGCCAATCAGTGTGAGAAGAGAAGACATGACAGAGCCAGGTGGAGCTGGAGAGC
 CAGCATTATGGAAAAACACCAGTTCTAGCATGGAGCCGACTGCACCCCTGGTGACTCCTCCACCCAAAGG
 AGGAGGCGACATGGCTGTGGTAGTGTCAAAGGAAGGTTCTGGGAGAAACCTAGTGATGACAGCTTTCAG
 AAGTCTGAAGCCAGGCCATCCCAGAGATGCCCATGTTTGAATCCCCAGTCTGACTTCAGTTTTCATG
 CTGATGAGTACCTAGAAGTCTACGTTTCTGCTTCTGAGCACCTAACCCTTCTGGATCCAGATCGTTGG
 CTCCCGCAGCCTGCAATTGGATAAGCTTGCAATGAGATGACCCAGCACTATGAGAATAGTGTGCCTGAA
 GACTTGACTGTGCATGTAGGAGACATTGTAGCAGCACCTTACCTACAAATGGTTCCTGGTATCGAGCCC
 GGGTCTCGGCACCTTGGAGAATGGGAACCTGGACCTCTATTTTGTGACTTTGGAGATAATGGAGATTG
 CCCACTGAAGGACCTCAGGGCTCTCAGGAGTGACTTCTAAGCCTTCCATTTCAAGCAATAGAATGTAGT
 CTGGCACGGATTGCTCCCTCAGGTGACCAGTGGGAAGAGGAAGCTTTGGATGAGTTTGATAGACTCACTC
 ATTGTGCTGACTGGAAGCCTCTGGTAGCCAAGATCTCTAGCTATGTCCAGACTGGGATCTCAACTTGGCC
 AAAGATCTACTTATGATACTAGCAATGGGAAGAACTTGATATTGGGCTAGAATTAGTACACAAAGGA
 TACGCAATTGAGCTTCTGAAGACATAGAAGAAAACAGAGCTGTCCAGACATGTTGAAGGACATGGCCA
 CAGAAACAGATGCCTCTCTCAGCACGTTGCTCACTGAGACAAAAGAGCTCTGGAGAGATAACACATAC
 CCTGTCCTGCCTCAGCTTATCAGAAGCTGCTTCCATGTCTGGTGTGATAACCTGAAGATGACTACTTA
 CTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC222658 representing NM_001083965
 Red=Cloning site Green=Tags(s)

MSTERTSWTSLSTIQKIALGLGIPASATVAYILYRRYRESREERLTFVGEDDIEIEMRVPQEAVKLIIGR
 QGANIKQLRKQTGARIDVDTEVDGERVLLISGFVQVCKAKAAIHQILTENTPVSEQLSVPQRSVGRRI
 GRGGETIRSICKASGAKITCDKESEGTLLL SRL IKISGTQKEVAAAKHLILEKVSEDEELRKRIHSAET
 RVPRKQPI SVRREDMTEPGGAGEPALWKNTSSSMEPTAPLVTPPKGGDMAVVVSKEGSWEKPSDDSFQ
 KSEAQAIPMPMFEIPSPDFSFADEYLEVYVSASEHPNHFVIQIVGSRSLQLDKLVNEMTQHYENSVPE
 DLTVHVGDIVAAPLPTNGSWYRVRVLTLENGLDLYFVDFDNGDCPLKDLRALRSDFLSLPFQAI ECS
 LARIAPSGDQWEEELDEFDRLTHCADWKPLVAKISSYVQTGI STWPKIYLYDTSNGKLDIGLELVHKG
 YAIELPEDIEENRAVPDMLKDMATETDASLSTLLTETKSSGEITHTLSCLSLSEAA SMSGDDNLEDDYL
 L

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8067_d03.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001083965

ORF Size: 1683 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001083965.2](#)

RefSeq Size: 2854 bp

RefSeq ORF: 1686 bp

Locus ID: 11022

UniProt ID: [Q9Y2W6](#)


Cytogenetics: 1q21.3

Protein Families: Transmembrane

MW: 61.9 kDa

Gene Summary: Participates in the primary piRNA biogenesis pathway and is required during spermatogenesis to repress transposable elements and prevent their mobilization, which is essential for the germline integrity. The piRNA metabolic process mediates the repression of transposable elements during meiosis by forming complexes composed of piRNAs and Piwi proteins and govern the methylation and subsequent repression of transposons. Required for the final steps of primary piRNA biogenesis by participating in the processing of 31-37 nt intermediates into mature piRNAs. May act in pi-bodies and piP-bodies by transferring piRNA precursors or intermediates to or between these granules.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC222658