

## Product datasheet for **RC222661L4V**

### TFIIS (TCEA2) (NM\_003195) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TFIIS (TCEA2) (NM_003195) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TFIIS
Synonyms:	TFIIS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003195
ORF Size:	897 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222661).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003195.4</a>
RefSeq Size:	1749 bp
RefSeq ORF:	900 bp
Locus ID:	6919
UniProt ID:	<a href="#">Q15560</a>
Cytogenetics:	20q13.33
Domains:	TFIIS, TFS2N, TFS2M
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 33.6 kDa

**Gene Summary:** The protein encoded by this gene is found in the nucleus, where it functions as an SII class transcription elongation factor. Elongation factors in this class are responsible for releasing RNA polymerase II ternary complexes from transcriptional arrest at template-encoded arresting sites. The encoded protein has been shown to interact with general transcription factor IIB, a basal transcription factor. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]