

## Product datasheet for **RC222663L2V**

### SP5 (NM\_001003845) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SP5 (NM_001003845) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SP5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001003845
ORF Size:	1194 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222663).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001003845.1</a>
RefSeq Size:	1941 bp
RefSeq ORF:	1197 bp
Locus ID:	389058
UniProt ID:	<a href="#">Q6BEB4</a>
Cytogenetics:	2q31.1
Protein Families:	Transcription Factors
MW:	41.8 kDa



[View online »](#)

**Gene Summary:**

Binds to GC boxes promoters elements. Probable transcriptional activator that has a role in the coordination of changes in transcription required to generate pattern in the developing embryo (By similarity).[UniProtKB/Swiss-Prot Function]