

Product datasheet for **RC222665**

HIF3 alpha (HIF3A) (NM_152794) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HIF3 alpha (HIF3A) (NM_152794) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	HIF3 alpha
Synonyms:	bHLHe17; HIF-3A; HIF3-alpha-1; IPAS; MOP7; PASD7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC222665 representing NM_152794
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGACTGGCAAGACCACAGGTCGACCACGGAGCTGCGCAAGGAAAAGTCCCGGGATGCGGCCCGCAGCC
GGCGCAGCCAGGAGACCGAGGTGCTGTACCAGCTGGCTCACACGCTGCCCTTCGCCCGCGGGCTCAGCGC
CCACCTGGACAAGGCCTCTATCATGCGCCTCACCATCAGCTACCTGCGCATGCACCGCCTCTGCGCCGCA
GGGGAGTGAACAGGTGGGAGCAGGGGAGAACCCTGGATGCCTGCTACCTGAAGGCCCTGGAGGGCT
TCGTCATGGTGCTCACCGCCGAGGGAGACATGGCTTACCTGTCGGAGAATGTCAGCAAACACCTGGGCCCT
CAGTCAGCTGGAGCTCATTGGACACAGCATCTTTGATTTTCATCCACCCCTGTGACCAAGAGGAGCTTCAG
GACGCCCTGACCCCCAGCAGACCCTGTCCAGGAGGAAGGTGGAGGCCCCACGGAGCGGTGCTTCTCCT
TGCGCATGAAGAGTACACTCACCAGCCGCGGGCGCACCCCTCAACCTCAAGGCGGCCACCTGGAAGGTGCT
GAACTGCTCTGGACATATGAGGGCTACAAGCCACCTGCGCAGACTTCTCCAGCTGGGAGCCCTGACTCA
GAGCCCCCGCTGCAGTGCCTGGTGTCTATCTGCGAAGCCATCCCCACCCAGGCAGCCTGGAGCCCCAC
TGGGCCGAGGGGCTTCCTCAGCCGCCACAGCCTGGACATGAAGTTCACCTACTGTGACGACAGGATTGC
AGAAGTGGCTGGCTATAGTCCCAGTACCTGATCGGCTGTTCCGCCTACGAGTACATCCACGCGCTGGAC
TCCGATGCGGTGACGAAGAGCATCCACACCTTGTGAGCAAGGGCCAGGCAGTAACAGGGCAGTATCGCT
TCTGGCCCGGAGTGGTGGCTACCTGTGGACCCAGACCAGGCCACAGTGGTGTGAGGGGACGGGGCC
CCAGTCGGAGAGTATCGTCTGTGTCCATTTTTAATCAGCCAGGTGGAAGAGACCGGAGTGGTGTGTCC
CTGGAGCAAACGGAGCAACACTCTGCAGACCCATTGAGCGGGCGCCCCCTCTCAGAAGGACACCCCTA
ACCTGGGGACAGCCTTGACACCCCTGGCCCCGGATCCTTGCCCTTCTGCACCCGCTTCCCTGAGCGA
GGCTGCCCTGGCCGTGACCCCCGCCGTTTCTGCAGCCCTGACCTCCGTGCGCTCCTGGGACCCATCCTG
GATGGGGCTTCAGTAGCAGCCACTCCCAGCACCCCGCTGGCCACACGGCACCCCCAAAGTCTCTTCGG
CTGATCTCCAGATGAACTACCTGTGGCACCAGAAATGTGCACAGACTTTCACCTCCGGGAAAGACAC
TGAGGCAGTGGAGACAGATTTAGATATAGCTCAGGATGCTGATGCTCTGGATTTGGAGATGCTGGCCCC
TACATCTCCATGGATGATGACTTCCAGCTCAACGCCAGCGAGCAGCTACCCAGGGCCTACCACAGACCTC
TGGGGGCTGTCCCCGGCCCCGTGCTCGGAGCTTCCATGGCTGTACCTCCAGCCCTTGAGCCCTCCCT
GCTACCCCGCTGGGGAGTGACCCCGGCTGAGCTGCTCCAGCCCTTCCAGAGGGGACCCCTCAGCATCC
TCTCCATGGCTGGGGCTCGGAAGAGGACCTGGCCCAGAGCTCAGAGGACGAGGACGAGGGAGTGGAGC
TGCTGGGAGTGGAGACTCCAAAAGTCCCCAGCCAGAACACGAAAACTTTCTGCTCTTCTCCTCAG
CCTGAGTTTCTTCTGACAGGAGGACCAGCCCCAGGGAGCCTGCAGGACCCAGCACCCACTCCTGAAC
CTGAATGAGCCCCCTGGGCTGGGCCCTCACTGCTCTCTCCGTACTCAGACGAGGACACTACCCAGCCCC
GGGGCCCCCTCCAGCCAAGGGCAGGCTCAGCCCAGGCTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222665 representing NM_152794
Red=Cloning site Green=Tags(s)


MDWQDHRSTTELKREKSRDAARRRSQETEVLYQLAHTLPFARGVSAHLKASIMRLTISYLRMHRCAA
 GEWNQVVGAGGEPLDACYLKALEGFVMVLTAEGDMAYLSENVSKHLGLSQLELIGHISIFDFIHPCDQEELQ
 DALTPQQTLSRRKVEAPTERCFSLRMKSTLTSRGRTLNLKAATWKVNLNCSGHRMAYKPPAQTSPAGSPDS
 EPPLQCLVLICEAIPHPGSLPPPLGRGAFLSRHSLDMKF TYCDDRIA EVAGYSPDDLIGCSAYEYIHALD
 SDAVSKSIHTLLSKGQAVTGQYRFLARSGYLWTQTQATVVSGGRGPQSEI VCVHFLISQVEETGVVLS
 LEQTEQHSRRPIQRGAPSQKDTNPNGDSLDTPGPRILAF LHPPSLSEAALAADPRRF CSPDLRRLGPIIL
 DGASVAATPSTPLATRHPQSPLSADLPDEL PVGTENVHRLFTSGKDTEAVETDL DIAQDADALDLEMLAP
 YISMDDDFQLNASEQLPRAYHRPLGAVPRPRARSFHGLSPPALEPSLLPRWGS DPRLSCSSPSRGPDSAS
 SPMAGARKRTLAQSSEDEDEGVELLGVRPPKRSPEHENFLLFPLSLFLLTGGPAPGSLQDPSTPLLN
 LNEPLGLGPSLLSPYSDEDTTQPGGPFQPRAGSAQAD

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8072_b06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_152794

ORF Size: 2001 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152794.4
RefSeq Size:	2717 bp
RefSeq ORF:	2004 bp
Locus ID:	64344
UniProt ID:	Q9Y2N7
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome, Transcription Factors
MW:	72.2 kDa
Gene Summary:	<p>The protein encoded by this gene is the alpha-3 subunit of one of several alpha/beta-subunit heterodimeric transcription factors that regulate many adaptive responses to low oxygen tension (hypoxia). The alpha-3 subunit lacks the transactivation domain found in factors containing either the alpha-1 or alpha-2 subunits. It is thought that factors containing the alpha-3 subunit are negative regulators of hypoxia-inducible gene expression. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2011]</p>

Product images:



Circular map for RC222665