

Product datasheet for RC222670

WNT8B (NM 003393) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: WNT8B (NM_003393) Human Tagged ORF Clone

Tag:Myc-DDKSymbol:WNT8B

Mammalian Cell Neomycin

Selection:

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

ORF Nucleotide >RC222670 representing NM_003393

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGTTTCTTTCAAAGCCTTCTGTGTACATCTGTCTTTTCACCTGTGTCCTCCAACTCAGCCACAGCTGGT CGGTGAACAATTTCCTGATGACTGGTCCAAAGGCTTACCTGATTTACTCCAGCAGTGTGGCAGCTGGTGC CCAGAGTGGTATTGAAGAATGCAAGTATCAGTTTGCCTGGGACCGCTGGAACTGCCCTGAGAGAGCCCTG CAGCTGTCCAGCCATGGTGGGCTTCGCAGTGCCAATCGGGAGACAGCATTTGTGCATGCCATCAGTTCTG CTGGAGTCATGTACACCCTGACTAGAAACTGCAGCCTTGGAGATTTTGATAACTGTGGCTGTGATGACTC CCGCAACGGGCAACTGGGGGGACAAGGCTGGCTGTGGGGAGGCTGCAGTGACAATGTGGGCTTCGGAGAG GCGATTTCCAAGCAGTTTGTCGATGCCCTGGAAACAGGACAGGATGCACGGGCAGCCATGAACCTGCACA ACAACGAGGCTGGCCGCAAGGCGGTGAAGGGCACCATGAAACGCACGTGCAAGTGCCACGGCGTGTCTGG CAGCTGCACCACGCAGACCTGTTGGCTGCAGCTGCCCGAGTTCCGCGAGGTGGGCGCGCACCTGAAGGAG TCGCCGACACCTTTCGCTCCATCTCTACCCGGGAGCTGGTGCACCTGGAGGACTCCCCGGACTACTGCCT GGAGAACAAAACGCTAGGGCTGCTGGGCACCGAAGGCCGAGAGTGCCTAAGGCGCGGGCCGGGCCCTGGGT CGCTGGGAACGCCGCAGCTGCCGCCGGCTCTGCGGGGACTGCGGGCTGGAGGAGCGCCGGGCCG AGACCGTGTCCAGCTGCAACTGCAAGTTCCACTGGTGCTGCGCAGTCCGCTGCGAGCAGTGCCGCCGGAG CCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Protein Sequence: >RC222670 representing NM_003393

Red=Cloning site Green=Tags(s)

MFLSKPSVYICLFTCVLQLSHSWSVNNFLMTGPKAYLIYSSSVAAGAQSGIEECKYQFAWDRWNCPERAL QLSSHGGLRSANRETAFVHAISSAGVMYTLTRNCSLGDFDNCGCDDSRNGQLGGQGWLWGGCSDNVGFGE AISKQFVDALETGQDARAAMNLHNNEAGRKAVKGTMKRTCKCHGVSGSCTTQTCWLQLPEFREVGAHLKE KYHAALKVDLLQGAGNSAAGRGAIADTFRSISTRELVHLEDSPDYCLENKTLGLLGTEGRECLRRGRALG RWERRSCRRLCGDCGLAVEERRAETVSSCNCKFHWCCAVRCEQCRRRVTKYFCSRAERPRGGAAHKPGRK P


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6049 g02.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_003393

ORF Size: 1053 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: <u>NM 003393.4</u>

 RefSeq Size:
 2117 bp

 RefSeq ORF:
 1056 bp

 Locus ID:
 7479

 UniProt ID:
 Q93098

 Cytogenetics:
 10q24.31

Domains: wnt

Protein Families: Cancer stem cells, ES Cell Differentiation/IPS, Secreted Protein, Stem cell relevant signaling -

Wnt Signaling pathway

Protein Pathways: Basal cell carcinoma, Hedgehog signaling pathway, Melanogenesis, Pathways in cancer, Wnt

signaling pathway

MW: 38.5 kDa


Gene Summary: The WNT gene family consists of structurally related genes which encode secreted signaling

proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 95%, 86% and 71% amino acid identity to the mouse, zebrafish and Xenopus Wnt8B proteins, respectively. The expression patterns of the human and mouse genes appear identical and are restricted to the developing brain. The chromosomal location of this gene to 10q24 suggests it as a candidate


gene for partial epilepsy. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC222670



Western blot validation of overexpression lysate (Cat# [LY401156]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC222670 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).