

Product datasheet for **RC222673A1V**

Human DCAF4L2 (NM_152418) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DCAF4L2 (NM_152418) AAV Particle
Tag: Myc-DDK
Symbol: DCAF4L2
Synonyms: WDR21C
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC222673 representing NM_152418
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGAGCAAAAGACCGCGACTGCTCGAGGAAGCAGACAAGCAGAAAAAGACAGTCAGAGTGGGACTCA
ATGCACCTTCCATGCTACGAAAGAACCAGCTAGGTTTCCTCAGATTCGCCAACTATTGCCGTATAGCTCG
CGAGCTGCGTGTAAAGCTGCATGCAGAGGAAAAAGGTCCAGATTCATAGCTGGGATCCCTCCTCTTTGGCA
AGCGACCGATTAAACCGCATACTGGCGAATACCAACACTGACCAGCTCTTCACAGTGAACCAAGTCGAAG
CTGGAGGCTCCAAGTACGGCATCATCACCATGCGAGGCCCTGACGACCCCTGAGCTCCGGGTATACCCGCA
CAAAACCCTCTACGTCCCTAATCGGAAGGTGAATTCTATGTGCTGGGCCCTCACTGAATCACTTGGATTCC
CACCTTCTGCTGTGCTTCGTGGGACTTGACAGATACTCCAAGCTGTGCCGTGCTGCTCCAGCGTCTGCTGT
TCATAGGTAGCTTCCAGGAATGCGTCGGCCTGGCATGCTTTGCAGTTTCCAGATCCCTGATGCCTGGTC
CTGTGCCTGGTCCCTGAGCATCCACGCGTATCACTCTTCAGTACAGGCTTGTCTCAGCAGGTCTGTTG
ACCAACGTGGTGACGGGACACCAGCAGTCATTTGGGACTAGCAGTGATGTCTGGCCCAGCAGTTTGCAA
TCATGACTCCTTTGCTGTTAATGGCTGTGCTCTGGGAGATCTTTGGCATTGATCTGCGCTGTGGAAA
TCAAGGCAGCGGGTGAAGGCCATTTGCCTGTCCCATGATTAGCAGTGACTTCTCTGCAATCTCCAA
GATGGCCAATTCCTGGTGTATCAGACATGACTGGAATATCAAGCTGTGGGACTTGAGGGCCACTAAAT
GTGTAACACAGTACGAAGGTATGTGAATAACTCCGCCTACCTACCCGTGCATGTGAACGAAGAAGG
AGTCGTGGCGGCCGTGGGCCAGGACTGTACACGAGAATCTGGAGCCTCCGTATGGCCACCTGCTCACA
ACCATACCCTCCCATACCCCGCTCGGAGAACGACATTCACAGTGTGGCCTTCTCTTCTCGCCTCGGGG
GCTTCCGAGGAGCACCAGGGCTGCTCATGGCTGTCCGGGAGGACCTTTATTGTTTCTCTACGGT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC222673 representing NM_152418 Red=Cloning site Green=Tags(s) MESKRPRLLLEEADKQKKTVRVGLNAPSMRLRKNQLGFLRFANYCRIARELRVSCMQRKQVQIHSWDPSSLA SDRFNRILANTNTDQLFTVNQVEAGGSKYGIITMRGLTPELRVYPHKTLVVPNRKVNSMCWASLNHLDS HLLLCFVGLADTPSCAVLLPASLFIGSFPGMRRPGMLCSFQIPDAWSCAWSLSIHAYHSFSTGLSQQVLL TNVVTGHQQSFGTSSDVLAQQFAIMTPLLFGCRSGEIFGIDLRCGNQGSQGWKAICLSHDSAVTSLQILQ DGQFLVSSDMTGTIKLWDLRATKCVTQYEGHVNNSAYLPVHVNEEEGVVAAGVQDCYTRIWSLRHGHLTT TIPSPYPASENDIPSAFSSRLGGFRGAPGLLMAVREDLYCFSYG SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152418
ORF Size:	1185 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_152418.2
RefSeq Size:	3320 bp
RefSeq ORF:	1188 bp
Locus ID:	138009
UniProt ID:	Q8NA75
Cytogenetics:	8q21.3
MW:	43.6 kDa