

Product datasheet for **RC222721L1V**

HOXA1 (NM_153620) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HOXA1 (NM_153620) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOXA1
Synonyms:	BSAS; HOX1; HOX1F
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_153620
ORF Size:	331 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222721).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153620.2
RefSeq Size:	2358 bp
RefSeq ORF:	414 bp
Locus ID:	3198
UniProt ID:	P49639
Cytogenetics:	7p15.2
Protein Families:	Druggable Genome, Transcription Factors
MW:	14.6 kDa



[View online »](#)

Gene Summary:

In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. The encoded protein may be involved in the placement of hindbrain segments in the proper location along the anterior-posterior axis during development. Two transcript variants encoding two different isoforms have been found for this gene, with only one of the isoforms containing the homeodomain region. [provided by RefSeq, Jul 2008]