

Product datasheet for **RC222765L3V**

ADAM23 (NM_003812) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAM23 (NM_003812) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAM23
Synonyms:	MDC-3; MDC3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003812
ORF Size:	2496 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222765).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003812.2
RefSeq Size:	3059 bp
RefSeq ORF:	2499 bp
Locus ID:	8745
UniProt ID:	O75077
Cytogenetics:	2q33.3
Domains:	Reprolysin, DISIN, Pep_M12B_propep, ACR
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane



[View online »](#)

MW: 91.7 kDa

Gene Summary: This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. It is reported that inactivation of this gene is associated with tumorigenesis in human cancers. [provided by RefSeq, May 2013]