

## Product datasheet for **RC222794L4V**

### ZNF645 (CBLL2) (NM\_152577) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF645 (CBLL2) (NM_152577) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CBLL2
Synonyms:	CT138; HAKAIL; ZNF645
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152577
ORF Size:	1275 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222794).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152577.1</a>
RefSeq Size:	1560 bp
RefSeq ORF:	1278 bp
Locus ID:	158506
UniProt ID:	<a href="#">Q8N7E2</a>
Cytogenetics:	Xp22.11
MW:	48.8 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a member of the zinc finger domain-containing protein family. This family member contains both a RING-type and a C2H2-type of zinc finger domain, and it may function as an E3 ubiquitin-protein ligase. Protein localization suggests a role in human sperm production and quality control. [provided by RefSeq, Aug 2011]