

## Product datasheet for **RC222824L3V**

### SLC25A3 (NM\_005888) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC25A3 (NM_005888) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC25A3
Synonyms:	OK/SW-cl.48; PHC; PTP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005888
ORF Size:	1086 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222824).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005888.2</a>
RefSeq Size:	1688 bp
RefSeq ORF:	1089 bp
Locus ID:	5250
UniProt ID:	<a href="#">Q00325</a>
Cytogenetics:	12q23.1
Domains:	mito_carr
Protein Families:	Transmembrane



[View online »](#)

**MW:** 40.09 kDa

**Gene Summary:** The protein encoded by this gene catalyzes the transport of phosphate into the mitochondrial matrix, either by proton cotransport or in exchange for hydroxyl ions. The protein contains three related segments arranged in tandem which are related to those found in other characterized members of the mitochondrial carrier family. Both the N-terminal and C-terminal regions of this protein protrude toward the cytosol. Multiple alternatively spliced transcript variants have been isolated. [provided by RefSeq, Jul 2008]