

Product datasheet for **RC222832A1V**

Human DAPL1 (NM_001017920) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DAPL1 (NM_001017920) AAV Particle
Tag: Myc-DDK
Symbol: DAPL1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC222832 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

**ATGGCAAATGAAGTGCAAGACCTGCTCTCCCCTCGAAAGGGGGACATCCTCCTGCAGTAAAAGCTGGAG
GAATGAGAATTTCCAAAAACAAGAAATGGCACCTTGGAAAGACATACCAAAAAACAGGATTCGAGAA
AACAAAGTGCCATTGCAAATGTTGCCAAAATACAGACACCGGATGCCCTGAATGACGCACTGGAGAAGCTC
AACTATAAATTTCCAGCAACAGTGCACATGGCACATCAAAAACCCACACCTGCTCTGGAAAAGGTTGTC
CACTGAAAAGGATCTACATTATTCAGCAGCCTCGAAAATGT**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222832 protein sequence
Red=Cloning site **Green**=Tags(s)

MANEVQDLLSPRKGHPPAVKAGMIRISKKQEIGTLERHTKKTGFEKTSIAIANVAKIQTPDALNDALEKL
NYKFPATVHMAHQKPTPALEKVVPLKRIYIIQQPRKC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001017920
ORF Size: 321 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_001017920.1](#)

RefSeq Size: 551 bp

RefSeq ORF: 324 bp

Locus ID: 92196

UniProt ID: [A0PJW8](#)

Cytogenetics: 2q24.1

MW: 11.9 kDa