

Product datasheet for **RC222850L3V**

C10orf46 (CACUL1) (NM_153810) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C10orf46 (CACUL1) (NM_153810) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C10orf46
Synonyms:	C10orf46; CAC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153810
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222850).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153810.4
RefSeq Size:	6758 bp
RefSeq ORF:	1110 bp
Locus ID:	143384
UniProt ID:	Q86Y37
Cytogenetics:	10q26.11
MW:	40.9 kDa


[View online »](#)

Gene Summary: Cell cycle associated protein capable of promoting cell proliferation through the activation of CDK2 at the G1/S phase transition.[UniProtKB/Swiss-Prot Function]