

## Product datasheet for **RC222859L1V**

### **HOPX (NM\_032495) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HOPX (NM_032495) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOPX
Synonyms:	CAMEO; HOD; HOP; LAGY; NECC1; OB1; SMAP31; TOTO
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_032495
ORF Size:	219 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222859).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032495.4</a>
RefSeq Size:	1552 bp
RefSeq ORF:	276 bp
Locus ID:	84525
UniProt ID:	<a href="#">Q9BPY8</a>
Cytogenetics:	4q12
Domains:	homeobox
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 8.3 kDa

**Gene Summary:** The protein encoded by this gene is a homeodomain protein that lacks certain conserved residues required for DNA binding. It was reported that choriocarcinoma cell lines and tissues failed to express this gene, which suggested the possible involvement of this gene in malignant conversion of placental trophoblasts. Studies in mice suggest that this protein may interact with serum response factor (SRF) and modulate SRF-dependent cardiac-specific gene expression and cardiac development. Multiple alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Feb 2009]