

## Product datasheet for **RC222876L3V**

### **SERCA3 (ATP2A3) (NM\_174956) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SERCA3 (ATP2A3) (NM_174956) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SERCA3
Synonyms:	SERCA3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_174956
ORF Size:	3087 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222876).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_174956.2</a> , <a href="#">NP_777616.1</a>
RefSeq Size:	4840 bp
RefSeq ORF:	3090 bp
Locus ID:	489
UniProt ID:	<a href="#">Q93084</a>
Cytogenetics:	17p13.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Alzheimer's disease, Calcium signaling pathway



[View online »](#)

**MW:** 112.5 kDa

**Gene Summary:** This gene encodes one of the SERCA Ca(2+)-ATPases, which are intracellular pumps located in the sarcoplasmic or endoplasmic reticula of muscle cells. This enzyme catalyzes the hydrolysis of ATP coupled with the translocation of calcium from the cytosol to the sarcoplasmic reticulum lumen, and is involved in calcium sequestration associated with muscular excitation and contraction. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]