

## Product datasheet for **RC222880**

### OR51A2 (NM\_001004748) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	OR51A2 (NM_001004748) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR51A2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC222880 representing NM_001004748 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCATTATCAACACATCATATGTTGAAATCACCCACCTTCTTCTTGGTTGGGATGCCAGGGCTAGAAT  
ATGCACACATCTGGATCTCTATCCCCATCTGCAGCATGTATCTTATTGCTATTCTAGGAAATGGCACCAT  
TCTTTTTATCATCAAGACAGAGCCCTCCTTGCATGGGCCCATGTACTATTTTCTTCCATGTTGGCTATG  
TCAGACTTGGGTTTGTCTTTATCATCTCTGCCACTGTGTTAAGCATCTTCTGTTCAATGCCCTGAAA  
CTTCTTCTAGTGCCTGCTTTGCCAGGAATCTTCATTCATGGATTCTCAGTACTGGAGTCTCAGTCCT  
CCTGATCATGTCAATTTGATAGATTCCCTAGCCATCCACAATCCTCTGAGATACACCTCAATCCTGACAACT  
GTCAGAGTTGCCCAAATAGGGATAGTATTCTCCTTTAAGAGCATGCTCCTGGTCTTCCCTTCCCTTTCA  
CTTTAAGAAGCTTGAGATATTGCAAGAAAAACCAATTATCCCATTCTACTGTCTCCACCAGGATGTCAT  
GAAGTTGGCCTGTTCTGACAACAGAATTGATGTTATCTATGGCTTTTTTGGAGCACTCTGCCTTATGGTA  
GACTTTATTCTATTGCTGTGTCTTACACCCTGATCCTCAAGACTGTACCGGAATTGCATCCAAAAAGG  
AGGAGCTAAGGCTCTCAATACTGTGTTTCACACATCTGTGCAGTGATCATCTTCTACCTGCCATCAT  
CAACCTGGCCGTTGTCCACCGCTTGGCCGGCATGTCTCTCCCCTATTAATGTTCTCATGGCAAATGTT  
CTCCTACTGTACCTCCGCTGATGAAACCAATTGTTTATTGTGTAATAACTAAACAGATTAGAGTGAGAG  
TTGTAGCAAATTTGTGCAATGGAAGATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC222880 representing NM\_001004748  
Red=Cloning site Green=Tags(s)

MSIINTSYVEITTFVVGMPGLEAYAHIWISIPICSMYLIAILGNGTILFIIKTEPSLHGPMYYFLSMLAM  
 SDLGLSLSSLPTVLSIFLFPAPETSSSACFAQEFFIHGFVLESSVLLIMSFDRFLAIHNPLRYTSILTT  
 VRVAQIGIVFSFKSMLLVLPFPFTLRSLRYCKKNQLSHSYCLHQDVMKLCASDNRIDVIYGFAGALCLMV  
 DFILIAVSYTLILKTVPGIASKKEELKALNTCVSHICAVIIFYLPIINLAVVHRFAGHVSPILINVLMANV  
 LLLVPPMLKPIVYCVKTKQIRVRVVAKLCQWKI


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8002\\_h03.zip](https://cdn.origene.com/chromatograms/mk8002_h03.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_001004748

**ORF Size:** 939 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001004748.1</a> , <a href="#">NP_001004748.1</a>
<b>RefSeq Size:</b>	942 bp
<b>RefSeq ORF:</b>	942 bp
<b>Locus ID:</b>	401667
<b>UniProt ID:</b>	<a href="#">Q8NGJ7</a>
<b>Cytogenetics:</b>	11p15.4
<b>Protein Families:</b>	Transmembrane
<b>Protein Pathways:</b>	Olfactory transduction
<b>MW:</b>	34.9 kDa
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC222880