

Product datasheet for **RC222896**

MAD1 (MAD1L1) (NM_001013837) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAD1 (MAD1L1) (NM_001013837) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MAD1
Synonyms:	MAD1; PIG9; TP53I9; TXBP181
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC222896 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAAGACCTGGGGAAAAACCATGGTTTTATCCACCCTGAGATCTTTGAACAACTTCATCTCTCAGC
 GTGTGGAGGGAGGCTCTGGACTGGATATTTCTACCTCGGCCCCAGGTTCTCTGCAGATGCAGTACCAGCA
 GAGCATGCAGCTGGAGGAAAGAGCAGAGCAGATCCGTTTGAAGTCCCACCTCATCCAGGTGGAGCGGGAG
 AAAATGCAGATGGAGCTGAGTCACAAGAGGGCTCGAGTGGAGCTGGAGAGAGCAGCCAGCACCAGTGCCA
 GGAACACGAGCGTGAGGTGACCCGAACAGGAGCTCCTGACGCGCATCCGGCAGCTTCAGGAGCGGGA
 GGCCGGGGCGGAGGAGAAGATGCAGGAGCAGCTGGAGCGCAACAGGCAGTGTGAGCAGAACTGGATGCT
 GCCAGCAAGAGGCTGCGTGAGAAAGAGGACAGTCTGGCCAGGCTGGCGAGACCATCAACGCACTGAAGG
 GGAGGATCTCGAACTGCAGTGGAGCGTATGGACCAGGAGATGCGGGTGAAGCGCCTGGAGTCGGAGAA
 GCAGGAGCTGCAGGAGCAGCTGGACCTGCAACACAAAAATGCCAGGAAGCCAATCAGAAAAATCCAGGAA
 CTCCAGGCCAGCCAAGAAGCAAGAGCAGACCACGAGCAGCAGATTAAGGATCTGGAGCAGAAGCTGTCCC
 TGCAAGAGCAGGATGCAGCGATTGTGAAGAACATGAAGTCTGAGCTGGTACGGCTCCCTAGGCTGGAACG
 GGAGCTGAAGCAGCTGCGGGAGGAGAGCGCACCTGCGGGAGATGAGAGAGACCAACGGGCTGCTCCAG
 GAAGAGCTGGAAGGGCTGCAGAGGAAGCTGGGGCGCCAGGAGAAGATGCAGGAGACGCTGGTTGGCTTG
 AGCTGGAGAACGAGAGGCTGCTGGCCAAGCTGCAAAGCTGGGAGAGACTGGACCAGACCATGGGCCTGAG
 CATCAGGACTCCAGAAGACCTTTCCAGATTCTGGTTGAGCTGCAGCAGAGGGAGCTTGCCTTGAAGGAC
 AAGAACAGCGCCGTCACCAGCAGCGCCCGGGGCTGGAGAAGGCCAGGCAGCAGCTGCAGGAGGAGCTCC
 GCGAGGTCAGCGCCAGCTGTTGGAGGAGAGGAAGAAGCCGAGACCCACGAGGCGCTGGCCCGGAGGCT
 CCAGAAAACGGGTCTGCTGCTACCAAGGAGCGGGACGGTATGCGGGCCATCCTGGGGTCTACGACAGC
 GAGCTGACCCCGCCGAGTACTACCCAGCTGACGCGGCGCATGCGGGAGGCTGAGGATATGGTGCAGA
 AGGTGCACAGCCACAGCGCCGAGATGGAGGCTCAGCTGTGCGAGGCCCTGGAGGAGCTGGGAGGCCAGAA
 ACAAAGAGCAGACATGCTGGAGATGGAGCTGAAGATGCTGAAGTCTCAGTCCAGCTCTGCCGAACAGAGC
 TTCCTGTTCTCCAGGGAGGAGGCGGACACGCTCAGGTTGAAGGTCGAGGAGCTGGAAGGCGAGCGGAGTC
 GGCTGGAGGAGGAAAAGAGGATGCTGGAGGCACAGCTGGAGCGGCGAGCTCTGCAGGGTACTATGACCA
 GAGCAGGACCAAAGTGTGCACATGAGCCTGAACCCACCAGTGTGCCAGGCAGCGCCTGCGCGAGGAC
 CACAGCCAGCTGCAGGCGGAGTGGGAGCGACTGCGCGGGCTCCTGCGCGCATGGAGAGAGGAGGCACCG
 TCCCAGCCGACCTTGAGGCTGCCGCCGAGTCTGCCATCGTCCAAGGAGGTGGCAGAGCTGAAGAAGCA
 GGTGGAGAGTGCCGAGCTGAAGAACCAGCGGCTCAAGGAGTTTTCCAGACCAAGATCCAGGAGTTCCGC
 AAGGCCTGCTACACGCTCACCGCTACCAGATCGACATCACACGGAGAACCAGTACCGGCTGACCTCGC
 TGTACGCCGAGCACCCAGGCGACTGCCTCATCTTCAAGGCCACCAGCCCCCTGGGTTCCAAGATGCAGCT
 ACTGGAGACAGAGTTCTCACACCCGTGGGCGAGCTCATCGAGGTGCACCTGCGGCGCCAGGACAGCATC
 CCTGCCTCTCAGCTCGCTCACCTCGAGCTCTCAGCCGCCAGACCGTGGCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222896 protein sequence
Red=Cloning site Green=Tags(s)

```
MEDLGENTMVLSTLRSLNNFISQRVEGGSGLDISTSAPGSLQMQYQQSMQLEERAEQIRSKSHLIQVERE
KMQMELSHKRARVELERAAS TSARNYEREVDRNQELL TRIRQLQEREAGAEKMQEQLERNRQCQQLDA
ASKRLREKEDSLAQAGETINALKGRISELQWSVMDQEMRVKRLSEKQELQEQLDLQHKKCQEANQKIQE
LQASQEARADHEQQIKDLEQKLSLQEQDAAIKNNMKSELVRLPRLERELKQLREESAHLREMRETNGLLQ
EELEGLQRKLGRQEKMQETLVGLELENERLLAKLQSWERLDQTMGLSIRTPEDLSRFVVELQQRELALKD
KNSAVTSSARGLEKARQQLQEELRQVSGQLLEERKKRETHEALARRLQKRVLLLTKERDGMRAILGSYDS
ELTPAEYSPQLTRRMREAEDMVQKVHSHSAEMEAQLSQALEELGGQKQRADMLEMELKMLKSQSSSAEQS
FLFSREEADTLRLKVEELEGEERSRLEEEKRMLEAQLERRALQGDYDQSRTKVLHMSLNPTSVARQRLRED
HSQLQAECERLRGLLRAMERGGTVPADLEAAAASLPSSKEVAELKKQVESAE LKNQRLKEVFQTKIQEFR
KACYTLTGYQIDITTENQYRLTSLYAEHPGDCLIFKATSPSGSKMQLLETEFSHTVGELIEVHLRRQDSI
PAFLSSLTLELFSRQTV
```



TRTRPLEQKLISEEDLANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6690_b07.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001013837


ORF Size: 2154 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013837.2</u>
RefSeq Size:	2714 bp
RefSeq ORF:	2157 bp
Locus ID:	8379
UniProt ID:	<u>Q9Y6D9</u>
Cytogenetics:	7p22.3
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle
MW:	83.1 kDa
Gene Summary:	MAD1L1 is a component of the mitotic spindle-assembly checkpoint that prevents the onset of anaphase until all chromosome are properly aligned at the metaphase plate. MAD1L1 functions as a homodimer and interacts with MAD2L1. MAD1L1 may play a role in cell cycle control and tumor suppression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]

Product images:



Circular map for RC222896



Western blot validation of overexpression lysate (Cat# [LY423038]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC222896 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).