

Product datasheet for **RC222899L3V**

SLC22A7 (NM_153320) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC22A7 (NM_153320) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC22A7
Synonyms:	hOAT11; NLT; OAT2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153320
ORF Size:	1644 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222899).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153320.2 , NP_696961.2
RefSeq Size:	2559 bp
RefSeq ORF:	1647 bp
Locus ID:	10864
UniProt ID:	Q9Y694
Cytogenetics:	6p21.1
Domains:	sugar_tr
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 60 kDa

Gene Summary: The protein encoded by this gene is involved in the sodium-independent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and appears to be localized to the basolateral membrane of the kidney. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]