

## Product datasheet for **RC222913L3V**

### GCOM1 (NM\_001018091) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GCOM1 (NM_001018091) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GCOM1
Synonyms:	gcom; Gcom2; GRINL1A; MYZAP; MYZAP-POLR2M
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001018091
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222913).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001018091.1</a>
RefSeq Size:	4459 bp
RefSeq ORF:	1338 bp
Locus ID:	145781
UniProt ID:	<a href="#">P0CAP1</a>
Cytogenetics:	15q21.3
Protein Families:	Druggable Genome
MW:	51.6 kDa



[View online »](#)

**Gene Summary:**

This locus represents naturally occurring readthrough transcription between the neighboring MYZAP (myocardial zonula adherens protein) and POLR2M (polymerase (RNA) II (DNA directed) polypeptide M) genes on chromosome 15. Alternative splicing results in multiple readthrough transcript variants. Readthrough variants may encode proteins that share sequence identity with the upstream gene product or with both the upstream and downstream gene products. Some readthrough transcript variants are also expected to be candidates for nonsense-mediated decay (NMD). [provided by RefSeq, Oct 2013]