

## Product datasheet for **RC222922L2V**

### SLC15A1 (NM\_005073) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC15A1 (NM_005073) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC15A1
Synonyms:	HPECT1; HPEPT1; PEPT1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005073
ORF Size:	2124 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222922).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005073.1</a>
RefSeq Size:	3142 bp
RefSeq ORF:	2127 bp
Locus ID:	6564
UniProt ID:	<a href="#">P46059</a>
Cytogenetics:	13q32.2-q32.3
Protein Families:	Transmembrane
MW:	78.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an intestinal hydrogen peptide cotransporter that is a member of the solute carrier family 15. The encoded protein is localized to the brush border membrane of the intestinal epithelium and mediates the uptake of di- and tripeptides from the lumen into the enterocytes. This protein plays an important role in the uptake and digestion of dietary proteins. This protein also facilitates the absorption of numerous peptidomimetic drugs. [provided by RefSeq, Apr 2010]