

Product datasheet for **RC222931L4V**

SKA3 (NM_145061) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SKA3 (NM_145061) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SKA3
Synonyms:	C13orf3; RAMA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145061
ORF Size:	1236 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222931).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145061.3 , NP_659498.2
RefSeq Size:	2907 bp
RefSeq ORF:	1239 bp
Locus ID:	221150
UniProt ID:	Q8IX90
Cytogenetics:	13q12.11
MW:	46.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a component of the spindle and kinetochore-associated protein complex that regulates microtubule attachment to the kinetochores during mitosis. The encoded protein localizes to the outer kinetechore and may be required for normal chromosome segregation and cell division. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]