

Product datasheet for **RC222939L3V**

ONECUT3 (NM_001080488) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ONECUT3 (NM_001080488) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ONECUT3
Synonyms:	OC3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001080488
ORF Size:	1482 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222939).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001080488.1
RefSeq Size:	1485 bp
RefSeq ORF:	1485 bp
Locus ID:	390874
UniProt ID:	O60422
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome
MW:	49.9 kDa



[View online »](#)

Gene Summary: Transcriptional activator. Binds the consensus DNA sequence 5'-DHWATTGAYTWWD-3' on a variety of gene promoters such as those of HNF3B and TTR (By similarity).[UniProtKB/Swiss-Prot Function]