

Product datasheet for **RC222952L3V**

DNA Ligase III (LIG3) (NM_002311) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNA Ligase III (LIG3) (NM_002311) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LIG3
Synonyms:	LIG2; LIG3alpha
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002311
ORF Size:	2847 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222952).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002311.4 , NP_002302.2
RefSeq Size:	3103 bp
RefSeq ORF:	2850 bp
Locus ID:	3980
UniProt ID:	P49916
Cytogenetics:	17q12
Domains:	DNA_ligase, zf-PARP
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Base excision repair

MW: 106 kDa

Gene Summary: This gene is a member of the DNA ligase family. Each member of this family encodes a protein that catalyzes the joining of DNA ends but they each have a distinct role in DNA metabolism. The protein encoded by this gene is involved in excision repair and is located in both the mitochondria and nucleus, with translation initiation from the upstream start codon allowing for transport to the mitochondria and translation initiation from a downstream start codon allowing for transport to the nucleus. Additionally, alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]