

Product datasheet for RC222962L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Nectin 3 (NECTIN3) (NM 015480) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Nectin 3 (NECTIN3) (NM 015480) Human Tagged ORF Clone Lentiviral Particle

Symbol: Nectin 3

Synonyms: CD113; CDW113; NECTIN-3; PPR3; PRR3; PVRL3; PVRR3

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_015480 **ORF Size:** 1647 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC222962).

Sequence:
OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 015480.1

RefSeq Size: 1650 bp
RefSeq ORF: 1650 bp
Locus ID: 25945
UniProt ID: Q9NQS3
Cytogenetics: 3q13.13

Domains: ig, IG

Protein Families: Druggable Genome, Transmembrane





Nectin 3 (NECTIN3) (NM_015480) Human Tagged ORF Clone Lentiviral Particle - RC222962L4V

Protein Pathways: Adherens junction, Cell adhesion molecules (CAMs)

MW: 60.8 kDa

Gene Summary: This gene encodes a member of the nectin family of proteins, which function as adhesion

molecules at adherens junctions. This family member interacts with other nectin-like proteins and with afadin, a filamentous actin-binding protein involved in the regulation of directional motility, cell proliferation and survival. This gene plays a role in ocular development involving the ciliary body. Mutations in this gene are believed to result in congenital ocular defects. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011]