

## Product datasheet for **RC222966L4V**

### USP33 (NM\_015017) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	USP33 (NM_015017) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP33
Synonyms:	VDU1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015017
ORF Size:	2826 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222966).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015017.3</a>
RefSeq Size:	4524 bp
RefSeq ORF:	2829 bp
Locus ID:	23032
UniProt ID:	<a href="#">Q8TEY7</a>
Cytogenetics:	1p31.1
Domains:	UCH, zf-UBP, DUSP
Protein Families:	Druggable Genome, Protease



[View online »](#)

**MW:** 106.5 kDa

**Gene Summary:** This gene encodes a deubiquinating enzyme important in a variety of processes, including Slit-dependent cell migration and beta-2 adrenergic receptor signaling. The protein is negatively regulated through ubiquitination by von Hippel-Lindau tumor protein (VHL). Alternative splicing results in multiple transcript variants and protein isoforms. [provided by RefSeq, Jun 2012]