

Product datasheet for **RC222981A1V**

Human TRAF4AF1 (KNSTRN) (NM_033286) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TRAF4AF1 (KNSTRN) (NM_033286) AAV Particle
Tag:	Myc-DDK
Symbol:	TRAF4AF1
Synonyms:	C15orf23; HSD11; SKAP; TRAF4AF1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC222981 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCTCCCGAAGCCCCGCCCTGGACAGAGTTTTCCGTACAACATGGCTGTCTACAGAGTGCATT
CCCACCCACTTCCGCCTAGCTACCGGAAGTTTCTATTTGAAACCCAGGAGGCCGACTTAGCCGGTGGCAC
GACAGTTGCTGCAGGGAATCTTTAAACGAGAGCGAGAAGGACTGCGGGCAGGACCGCGGGCTCCTGGG
GTTGAGCCGTGCCTCCTCGTTACGATGACCAAGTGTGGTTAAGACAGTGTATAGCCTGCAGCCCTCCTCTG
CGTGAGCGCGCGCCAGCCGGCAGACACAACTCGGGCCACTTCTAAGAGTCTCTTACCTGTTAGGTC
CAAAGAAGTCGATGTTTCCAAACAGCTTCATTAGGAGGTCAGAGAATGATGTTACAAAAATCACCAAA
CTGAGACGAGAGAATGGGCAATGAAAGCTACTGACACTGCCACCAGAAGGAATGTAGAAAAGGCTACA
AACCCTGAGTAAGCAAAAATCAGAGGAAGAGCTCAAGGACAAGAACCAGCTGTTAGAAGCCGTCAACAA
GCAGTTGCACCAGAAGTTGACTGAAACTCAGGGAGAGCTGAAGGACCTGACCCAGAAGGTAGAGCTGCTG
GAGAAGTTTCGGGACAACTGTTTGGCAATTTTGGAGAGCAAGGCGCTTGATCCAGCTTTAGGCGGTGAGA
CCCTGGCATCACGACAAGAATCCACTACTGATCACATGGACTCTATGTTGCTGTTAGAACTTTGAAGA
GGAGCTGAAGCTTTTAAACGAAACAGCCAAAAAGCAGATGGAGGAGTTACAGGCCTTAAAGGTAAGCTG
GAGATGAAAGAGGAAAGAGTCCGATTCTAGAACAGCAACCTTATGTAACAATCAAGTAAATGATTTAA
CAACAGCCCTTAAGGAAATGGAGCAGCTATTAGAAATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC222981 protein sequence Red=Cloning site Green=Tags(s) MAAPEAPPLDRVFRITWLSTECDSHPLPPSYRKFLFETQEADLAGGTTVAAGNLLNESEKDCGQDRRAPG VQPCLLVTMTSVVKTIVSLQPSSALSGGQPADTQTRATSKSLLPVRSEVDVSKQLHSGGPENDVTKITK LRRENGQMKATDTATRRNVKGYKPLSKQKSEELKDKNQLLEAVNKQLHQKL TETQGELKDLTQKVLL EKFRDNCLAILESKGLDPALGGETLASRQESTTDHMSMLLLETLEELKLFNETAKKQMEELQALKVKL EMKEERVRFLEQQTLCNNQVNDLTTALKEMEQLLEM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033286
ORF Size:	948 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033286.2, NP_150628.2</u>
RefSeq Size:	1763 bp
RefSeq ORF:	951 bp
Locus ID:	90417
UniProt ID:	<u>Q9Y448</u>
Cytogenetics:	15q15.1
MW:	35.4 kDa