

Product datasheet for **RC222985L2V**

SPTLC1 (NM_006415) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SPTLC1 (NM_006415) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPTLC1
Synonyms:	HSAN1; HSN1; LBC1; LCB1; SPT1; SPTI
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006415
ORF Size:	1419 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222985).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006415.2
RefSeq Size:	2780 bp
RefSeq ORF:	1422 bp
Locus ID:	10558
UniProt ID:	O15269
Cytogenetics:	9q22.31
Domains:	aminotran_1_2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: Metabolic pathways, Sphingolipid metabolism

MW: 52.6 kDa

Gene Summary: This gene encodes a member of the class-II pyridoxal-phosphate-dependent aminotransferase family. The encoded protein is the long chain base subunit 1 of serine palmitoyltransferase. Serine palmitoyltransferase converts L-serine and palmitoyl-CoA to 3-oxosphinganine with pyridoxal 5'-phosphate and is the key enzyme in sphingolipid biosynthesis. Mutations in this gene were identified in patients with hereditary sensory neuropathy type 1. Alternatively spliced variants encoding different isoforms have been identified. Pseudogenes of this gene have been defined on chromosomes 1, 6, 10, and 13. [provided by RefSeq, Jul 2013]