

Product datasheet for **RC223023**

OR5B17 (NM_001005489) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OR5B17 (NM_001005489) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: OR5B17
Synonyms: OR5B20P; OR11-237
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC223023 representing NM_001005489
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGAATAATACAGAGGTGAGTGAATTCATCCTGCTTGGTCTAACCAATGCCCCAGAACTACAGGTTC
CCCTCTTTATCATGTTTACCCTCATCTACCTCATCACTCTGACTGGAACCTGGGGATGATCATATTAAT
CCTGCTGGACTCTCATCTCCACTCCCATGTACTTTTTCTCAGTAACCTGTCTTGCAGGCATTGGT
TACTCCTCAGCTGCTCACTCCAAAGGTTTTAACTGGGTTGCTTATAGAAGACAAAGCCATCTCTACAGTG
CCTGTGCTGCTCAGATGTTCTTTGTGAGTCTTTGCCACTGTGAAAAATTACCTCTTGCCTCAATGGC
CTATGACCGCTACGCAGCAGTGTGTAACCCCTACATTATACCACCACCATGACAACACGTGTGTGTGCT
TGCTGGCTATAGGCTGTTATGTCATTGGTTTTCTGAATGCTTCTATCCAAATTGGAGATACATTTCCGC
TCTCTTTCTGCATGTCCAATGTGATTCACTTTTTCTGTGACAAACCAGCAGTCATTACTCTGACCTG
CTCTGAGAAACACATTAGTGAGTTGATTCTTGTCTTATCAAGTTTTAATGTCTTTTTGCACTTCTT
GTTACCTTGATTTCTATCTGTTTCATATTGATCACCATTCTTAAGAGGCACACAGGTAAGGGATACCAGA
AGCCTTTATCTACCTGTGGTTCTCACCTATTGCCATTTCTTATTTTATATAACTGTATCATCATGTA
CATACGACCAAGTCCAGTCATTCCATGGACACAGACAAAATTGCATCTGTGTTCTACACTATGATCATC
CCATGCTCAGTCCTATAGTCTATACCTGAGGAACAAAGACGTGAAGAATGCATTATGAAGTTGTTG
AGAAGGCAAAATATTCTCTAGATTCACTCTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC223023 representing NM_001005489
Red=Cloning site Green=Tags(s)


MENNTEVSEFILLGLTNAPELQVPLFIMFTLIYLITLTGNLGMIIILLDSLHPTMYFFLSNLSLAGIG
 YSSAVTPKVL TGLLIEDKAISSACAAQMFFCAV FAVTENYLLSSMAYDRYAAVCNPLHYTTMTTRVCA
 CLAIGCYVIGFLNASIQIGDTFRLSFCMSNVIIHFFCDKPAVITLTCSEKHISELILVLISSFNVFFALL
 VTLISYLFILITILKRHTGKG YQKPLSTCGSHLIAIFLFYITVIIMYIRPSSSHSMDTKIASVFYTMII
 PMLSPIVYTLRNKDVKN AFMKVVEKAKYSLDSVF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8004_a02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001005489

ORF Size: 942 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001005489.2](#)

RefSeq Size: 945 bp

RefSeq ORF: 945 bp

Locus ID: 219965

UniProt ID: [Q8NGF7](#)

Cytogenetics: 11q12.1


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 34.9 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC223023