

Product datasheet for **RC223026A1V**

Human OR5T1 (NM_001004745) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: OR5T1

Synonyms: OR5TIP; OR11-179

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223026 representing NM_001004745
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCAGGGTTGCCTTCAGATATGGATCTATACAAGCTTCAATTAACAATTTTACTGAAGTCACCATGT
TTATATTAATAAGCTTCACAGAAGAATTTGATGTGCAAGTCTTCCTATTTTATTATTTTAGCAATCTA
TCTATTCACCTCTAATAGGCAATTTAGGGCTGGTTGTACCGATCATTGGGGATTCTGGCTTCACAGCCCA
ATGTACTATTTTCTTGGTGTTCATTCTTGGATGTCTGCTATTCTACAGTTGTCACTCCAAAATGT
TGGTCAATTTCTGGCAAAAAATAATCTATTTCAATTTCTTGGATGTGCAACACAGATGTTTCTTGCTTG
TACTTTTGAACCAACAGAATGCTTTCTTGGCTGCAATGGCTTATGATCGCTATGTAGCCATCTACAAC
CCTCTCCTGTATTCAAGTGAAGCATGTACCCAGAGTCTATGTGCCACTCATCACTGCTTCCTATGTTGCTA
GCATTTTACATGCTACTATACATACAGTGGCTACATTTAGCCTGTCCTTCTGTGGATCCAATGAAATTAG
GCATGTCTTTTGAATATGCCTCCTCTGCTTGCTATTTCTTGTCTGACACTCACGTAATCCAGCTTCTA
TTCTTCTACTTTGTGGCTCTATTGAGATAGTCACTATCCTGATTGTCTGATCTCCTATGGTTTATTCT
TGTTGGCCATTCTGAAGATGCAGTCTGCTGAAGGGAGGAGAAAAGTCTTCTCTACATGTGGAGCTCACCT
AACTGGAGTGACAATTTATCATGGGACAATCCTCTTCAATGTATGTGAGACCAAGTTCAGCTACACTTCG
GACAATGACATGATAGTGTCAATTTTATACCATGTGATCCCATGCTGAATCCCATCATCTACAGTT
TGCGGAACAAAGATGTAAGGAGGCAATCAAAAGATTGCTTGTGAGAAATTGGTTCATAAATAAGTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC223026 representing NM_001004745 Red=Cloning site Green=Tags(s) MSGLPDMDLYKLQLNNFTEVTMFILISFTEEDVQVFLFLFLAIYLFTLIGNLGLVVPPIIGDFWLHSP MYYFLGVLSFLDVCYSTVVTPKMLVNFLAKNKSISFLGCATQMFLACTFGTTECFLLAAMAYDRYVAIYN PLLYSVMSPRVYVPLITASYVASILHATIHTVATFSLSFCGSNEIRHVCNMPLLAISCSDTHVIQLL FFYFVGSIEIVTILIVLISYGFILLAILKMQSAEGRRKVFSTCGAHLTGVTIYHGTILFMYVRPSSSYTS DNDMIVSIFYTIVIPMLNPIIYSLRNKDVKEAIKRLLRNWFINKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001004745
ORF Size:	978 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001004745.1, NP_001004745.1</u>
RefSeq Size:	981 bp
RefSeq ORF:	981 bp
Locus ID:	390155
UniProt ID:	<u>Q8NG75</u>
Cytogenetics:	11q12.1
MW:	36.7 kDa