

## Product datasheet for **RC223118L3V**

### WWP2 (NM\_007014) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	WWP2 (NM_007014) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WWP2
Synonyms:	AIP2; WWp2-like
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007014
ORF Size:	2610 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223118).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007014.3</a>
RefSeq Size:	4527 bp
RefSeq ORF:	2613 bp
Locus ID:	11060
UniProt ID:	<a href="#">O00308</a>
Cytogenetics:	16q22.1
Domains:	C2, HECT, WW
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Ubiquitin mediated proteolysis

**MW:** 98.7 kDa

**Gene Summary:** This gene encodes a member of the Nedd4 family of E3 ligases, which play an important role in protein ubiquitination. The encoded protein contains four WW domains and may play a role in multiple processes including chondrogenesis and the regulation of oncogenic signaling pathways via interactions with Smad proteins and the tumor suppressor PTEN. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 10. [provided by RefSeq, Jul 2012]