

## Product datasheet for RC223133L4

### UBE2V1 (TMEM189-UBE2V1) (NM\_199203) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	UBE2V1 (TMEM189-UBE2V1) (NM_199203) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	UBE2V1
Synonyms:	CROC-1B; CROC1B; KUA-UEV; TMEM189-UBE2V1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223133).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_199203
ORF Size:	1110 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_199203.1</a> , <a href="#">NP_954673.1</a>
<b>RefSeq Size:</b>	2852 bp
<b>RefSeq ORF:</b>	1113 bp
<b>Locus ID:</b>	387522
<b>UniProt ID:</b>	<a href="#">Q13404</a>
<b>Cytogenetics:</b>	20q13.13
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	42 kDa
<b>Gene Summary:</b>	The TMEM189-UEV mRNA is an infrequent but naturally occurring read-through transcript of the neighboring TMEM189 and UBE2V1 genes. Ubiquitin-conjugating E2 enzyme variant proteins constitute a distinct subfamily within the E2 protein family. They have sequence similarity to other ubiquitin-conjugating enzymes but lack the conserved cysteine residue that is critical for the catalytic activity of E2s. The protein produced by this transcript has UEV1 B domains but the protein is localized to the cytoplasm rather than to the nucleus. The significance of this read-through mRNA and the function of its protein product has not yet been determined. [provided by RefSeq, Oct 2010]

Product images:



Circular map for RC223133L4