

Product datasheet for **RC223183L4V**

RFP2 (TRIM13) (NM_005798) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RFP2 (TRIM13) (NM_005798) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RFP2
Synonyms:	CAR; DLEU5; LEU5; RFP2; RNF77
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005798
ORF Size:	1221 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223183).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005798.3
RefSeq Size:	6776 bp
RefSeq ORF:	1224 bp
Locus ID:	10206
UniProt ID:	O60858
Cytogenetics:	13q14.2
Domains:	zf-B_box, RING
Protein Families:	Transmembrane



[View online »](#)

MW: 47 kDa

Gene Summary: This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This gene is located on chromosome 13 within the minimal deletion region for B-cell chronic lymphocytic leukemia. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008]