

Product datasheet for **RC223195**

DBNDD2 (NM_001048222) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DBNDD2 (NM_001048222) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: DBNDD2
Synonyms: C20orf35; CK1BP; HSMNP1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC223195 representing NM_001048222
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

**ATGGACCCAAATCCTCGGGCCGCCCTGGAGCGCCAGCAGCTCCGCCTTCGGGAGCGGCAAAAATTCTTCG
AGGACATTTTACAGCCAGAGACAGAGTTTGTCTTCTCTGTCCCATCTGCATCTCGAGTCGAGAGACC
CCCCATAGGTAGTATCTCATCCATGGAAGTGAATGTGGACACACTGGAGCAAGTAGAACTTATTGACCTT
GGGACCCGGATGCAGCAGATGTGTTCTTGCCCTTGCGAAGATCCTCCACCAACCCCCAGTCGTCTGGTA
TGCCCTCTGCTTTGGGACTTCAGTGCCAGTCAGCCAGAGCCGGATGTCAGGCTC**

**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**

Protein Sequence: >RC223195 representing NM_001048222
Red=Cloning site Green=Tags(s)

MDPNPRAALERQQLRLRERQKFFEDILQPETEFVPLSHLHLESQRPPIGSISSEMEVNDTLEQVELIDL
GDPDAADVFLPCEDPPPTPQSSGMPLCFGDFSASQPEPDVRL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001048222

ORF Size: 336 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001048222.3</u>
RefSeq Size:	1263 bp
RefSeq ORF:	339 bp
Locus ID:	55861
UniProt ID:	<u>Q9BQY9</u>
Cytogenetics:	20q13.12
MW:	12.6 kDa
Gene Summary:	May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation (in vitro). [UniProtKB/Swiss-Prot Function]