

Product datasheet for **RC223231**

C16orf13 (METTL26) (NM_001040160) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C16orf13 (METTL26) (NM_001040160) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C16orf13
Synonyms: C16orf13; JFP2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC223231 representing NM_001040160
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGTGGCGGCGCCGCGGAGCGGAACAAGGATCCCATCTTGCACGTGCTGCGGCAGTACCTGGATC
CGGCCAGCGTGGCGTCCGCGTCTCGAGGTGGCCTCGGGCTCCGGCCAGCACGCAGCGCACTTCGCGCG
GGCCTTCCCCCTGGCCGAGTGGCAGCCGTCGGACGTGGACCAGCGCTGCCTGGACAGCATCGCGGCCACC
ACGCAAGCCCAGGGCCTGACCAACGTGAAGCCCCGCTACACCTGGACGTGACGTGGGGCTGGGAGCACT
GGGGCGGGATCTGCCACAGTCGCTGGACCTGTTGCTCTGCATCAACATGGCCCATGTCAGCCCCCTGCG
CTGCACGGAGGGGCTTTCAGAGCAGCAGGACACCTGCTCAAACCCAGGGCCCTGCTCATCACCTACGGG
GAACCCAGAATGGGGCTTCGGGACACAGCCCTCCTGGAGGACCTGGGAAAGGCCAGTGGCTGCTCCTG
GAGAGGATGGTGGACATGCCAGCCAACAACAAATGCCTGATCTTCGGAAAAACTAAGCCCTCCTTCAC
CCCCGCACACCTGCATCCCTGCCGAGGCTCTGTGAGGCACGAACCTGCCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223231 representing NM_001040160
Red=Cloning site Green=Tags(s)


MLVAAAERNKDPILHVLRLQYLDPAQRGVRVLEVASGSGQHAHFARAFPLAEWQPSDQVDRCLDSIAAT
TQAQGLTNVKAPLHLDVTWGWGILPQSLDLLLCINMAHVSPLRCTEGLFRAAGHLLKPRALLITYG
EPRMGASGHSPGGPGKGQWPAPGEDGGHASQQQMPDLPEKLSPSFTPAHLHPCRRRLCEARTLPP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001040160

ORF Size: 615 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001040160.3](#)

RefSeq Size: 786 bp

RefSeq ORF: 618 bp


Locus ID: 84326

UniProt ID: [Q96S19](#)

Cytogenetics: 16p13.3

MW: 22 kDa

Product images:



Circular map for RC223231