

## Product datasheet for **RC223261L3V**

### **KIR2.3 (KCNJ4) (NM\_152868) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	KIR2.3 (KCNJ4) (NM_152868) Human Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	KIR2.3
<b>Synonyms:</b>	HIR; HIRK2; HRK1; IRK-3; IRK3; Kir2.3
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
<b>Tag:</b>	Myc-DDK
<b>ACCN:</b>	NM_152868
<b>ORF Size:</b>	1335 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(RC223261).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_152868.1</a>
<b>RefSeq Size:</b>	2063 bp
<b>RefSeq ORF:</b>	1338 bp
<b>Locus ID:</b>	3761
<b>UniProt ID:</b>	<a href="#">P48050</a>
<b>Cytogenetics:</b>	22q13.1
<b>Protein Families:</b>	Druggable Genome, Ion Channels: Potassium, Transmembrane
<b>MW:</b>	49.3 kDa



[View online »](#)

**Gene Summary:**

Several different potassium channels are known to be involved with electrical signaling in the nervous system. One class is activated by depolarization whereas a second class is not. The latter are referred to as inwardly rectifying K<sup>+</sup> channels, and they have a greater tendency to allow potassium to flow into the cell rather than out of it. This asymmetry in potassium ion conductance plays a key role in the excitability of muscle cells and neurons. The protein encoded by this gene is an integral membrane protein and member of the inward rectifier potassium channel family. The encoded protein has a small unitary conductance compared to other members of this protein family. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]