

Product datasheet for **RC223278L3V**

DMRT2 (NM_181872) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DMRT2 (NM_181872) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DMRT2
Synonyms:	DSXL-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181872
ORF Size:	1683 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223278).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181872.1 , NP_870987.1
RefSeq Size:	1969 bp
RefSeq ORF:	1686 bp
Locus ID:	10655
UniProt ID:	Q9Y5R5
Cytogenetics:	9p24.3
Protein Families:	Transcription Factors, Transmembrane
MW:	61.6 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the DMRT gene family, sharing a DM DNA-binding domain with *Drosophila* 'doublesex' (dsx) and *C. elegans* mab3, genes involved in sex determination in these organisms. Also, this gene is located in a region of the human genome (chromosome 9p24.3) associated with gonadal dysgenesis and XY sex reversal. Hence this gene is one of the candidates for sex-determining gene(s) on chr 9. [provided by RefSeq, Apr 2010]