

Product datasheet for **RC223291L3V**

CCR5 (NM_000579) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CCR5 (NM_000579) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CCR5 |
| Synonyms: | CC-CKR-5; CCCR5; CCR-5; CD195; CKR-5; CKR5; CMKBR5; IDDM22 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_000579 |
| ORF Size: | 1056 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC223291). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_000579.3 , NP_000570.1 |
| RefSeq Size: | 3655 bp |
| RefSeq ORF: | 1059 bp |
| Locus ID: | 1234 |
| UniProt ID: | P51681 |
| Cytogenetics: | 3p21.31 |
| Protein Families: | Druggable Genome, ES Cell Differentiation/IPS, GPCR, Transmembrane |
| Protein Pathways: | Chemokine signaling pathway, Cytokine-cytokine receptor interaction, Endocytosis |




[View online »](#)


MW: 40.3 kDa

Gene Summary: This gene encodes a member of the beta chemokine receptor family, which is predicted to be a seven transmembrane protein similar to G protein-coupled receptors. This protein is expressed by T cells and macrophages, and is known to be an important co-receptor for macrophage-tropic virus, including HIV, to enter host cells. Defective alleles of this gene have been associated with the HIV infection resistance. The ligands of this receptor include monocyte chemoattractant protein 2 (MCP-2), macrophage inflammatory protein 1 alpha (MIP-1 alpha), macrophage inflammatory protein 1 beta (MIP-1 beta) and regulated on activation normal T expressed and secreted protein (RANTES). Expression of this gene was also detected in a promyeloblastic cell line, suggesting that this protein may play a role in granulocyte lineage proliferation and differentiation. This gene is located at the chemokine receptor gene cluster region. An allelic polymorphism in this gene results in both functional and non-functional alleles; the reference genome represents the functional allele. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2015]

Product images:



[RC223291L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC223291L3V particle to overexpress human CCR5-Myc-DDK fusion protein.



[RC223291L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC223291L3V particle to overexpress human CCR5-Myc-DDK fusion protein.