

Product datasheet for **RC223293**

MRPL55 (NM_181456) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MRPL55 (NM_181456) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: MRPL55
Synonyms: AAVG5835; L55nt; MRP-L55; PRO19675
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC223293 representing NM_181456
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCGGTGGGCAGCCTGCTTGGCCGGCTGAGGCAGAGCACCGTGAAGGCCACCGACCTGCACTCC
GCCGCTGCACACATCCTCCTGGCGAGCTGACAGCAGCAGGGCCCTCACTCACTCGTGTGCACCGCCAGGC
TTATGCACGACTCTACCCCGTGTGCTGGTGAAGCAGGATGGCTCCACCATCCACATCCGCTACAGGGAG
CCACGGCGCATGCTGGCGATGCCATAGATCTGGACACCCTGTCTCTGAGGAGCGCCGGCCAGGCTGC
GGAAGCGTGAGGCTCAGCTCCAGTCGAGGAAGGATACGAGCAGGAGCTCAGTGACTTGCATGTGGA
GCGCTACCGACAGTTCTGGACCAGGACCAAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223293 representing NM_181456
Red=Cloning site Green=Tags(s)
MAAVGSLLGRLRQSTVKATGPALRRLHTSSWRADSSRASLTRVHRQAYARLYPVLLVKQDGGSTIHIRYRE
PRRMLAMPIDLDTLSPEERRARLRKREAQLQSRKEYEQELSDDLHVERYRQFWTRTKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_181456.2 , NP_852121.1
RefSeq Size:	745 bp
RefSeq ORF:	387 bp
Locus ID:	128308
UniProt ID:	Q7Z7F7
Cytogenetics:	1q42.13
MW:	15.1 kDa
Gene Summary:	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Multiple transcript variants encoding two different isoforms were identified through sequence analysis. [provided by RefSeq, Jul 2008]