

Product datasheet for **RC223306L4V**

Keratinocyte differentiation associated protein (KRTDAP) (NM_207392) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Keratinocyte differentiation associated protein (KRTDAP) (NM_207392) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Keratinocyte differentiation associated protein
Synonyms:	KDAP; UNQ467
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207392
ORF Size:	297 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223306).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207392.1
RefSeq Size:	500 bp
RefSeq ORF:	300 bp
Locus ID:	388533
UniProt ID:	P60985
Cytogenetics:	19q13.12
MW:	11 kDa


[View online »](#)

Gene Summary:

This gene encodes a protein which may function in the regulation of keratinocyte differentiation and maintenance of stratified epithelia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]