

Product datasheet for **RC223315A1V**

Human **C11orf17 (AKIP1) (NM_020642) AAV Particle**

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: C11orf17

Synonyms: BCA3; C11orf17

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223315 representing NM_020642
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGACAAC**TGTTGGCGGCCG**CAGCGCTGAATGGGGTGGACCGACGTTCCCTGCAGCGTTCAGCAAGGC
TGGCTCTAGAAGTCTGGAGAGGGCCAAGAGGAGGGCGGTGGACTGGCATGCCCTGGAGCGTCCCAAAGG
CTGCATGGGGGTCTTGGCCGGAGGCGCCCCACCTAGAGAAACAGCCGGCAGCCGGCCCGCAGCGCGTT
CTCCCGGGAGAGAGAGAAGAGAGACCCCCAACCCCTTAGTGCTTCTTCAACAATGGCTGAATTCATGG
ACTATACTTCAAGTCAGTGTGGGAAATATTATTCATCTGTGCCAGAGGAAGGAGGGGCAACCCATGTCTA
TCGTTATCACAGAGGCGAGTCGAAGCTGCACATGTGCTTGGACATAGGGAATGGTCAGAGAAAAGACAGA
AAAAAGACATCCCTTGGTCTGGAGGCGAGCTATCAAAATACAGAGCATGCTCCAGAGGCATCCAGCGCTG
CTGAGAACATCTCTAAGGACCTCTACATAGAAGTATATCCAGGGACCTATTCTGTCACTGTGGGCTCAAA
TGACTTAACCAAGAAGACTCATGTGGTAGCAGTTGATTCTGGACAAAGCGTGGACCTGGTCTTCCCTGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAG**GTTTAA**

Protein Sequence: >RC223315 representing NM_020642
Red=Cloning site Green=Tags(s)

MDNCLAAAAALNGVDRSLQRSARLALEVLERAKRRAVDWHLERPKGCMGVLAAREAPHLEKQPAAGPQRV
LPGEREERPTLSASFRTMAEFMDYTSSQCGKYSSVPEEGGATHVYRYHRGESKLHMCLDIGNQQRKDR
KKTSLGPGGSYQISEHAPEASQPAENISKDLYIEVYPGTYSVTVGSNDLTKKTHVVAVDSGQSVDLVFPV

TRTRPLE**QKL**I**SEED**LAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_020642
ORF Size:	630 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_020642.2
RefSeq Size:	1332 bp
RefSeq ORF:	633 bp
Locus ID:	56672
UniProt ID:	Q9NQ31
Cytogenetics:	11p15.4
MW:	22.9 kDa